

Product datasheet for **MR201489L3V**

Tmem208 (NM_025486) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem208 (NM_025486) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem208
Synonyms:	1700006C06Rik; 2610030K20Rik; Hspc171
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025486
ORF Size:	522 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201489).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025486.1 , NP_079762.1
RefSeq Size:	791 bp
RefSeq ORF:	522 bp
Locus ID:	66320
UniProt ID:	Q9CR96
Cytogenetics:	8 D3
Gene Summary:	May function as a negative regulator of endoplasmic reticulum-stress induced autophagy. [UniProtKB/Swiss-Prot Function]



[View online »](#)