

## Product datasheet for MR201500

### Wfdc2 (BC099427) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Wfdc2 (BC099427) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Wfdc2  
**Synonyms:** HE4, WAP5  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201500 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCTGCCTGTCGCTCTGCTTGTGGCCGCTGGCCTCCTACTAGGGTTGCTACTGTTACCCCCATCT  
CAGCCACAGGCACCGATGCAGAGAAACCCGGCAGTGCCCCAGCTTGAACCAATTACGGACTGTGTGT  
GGAGTGCACTTTGGACAAGGACTGTGCGGACAACCGCAAGTGCTGCCAGGCGGGTGCAGCTCTGTCTGC  
TCCAAGCCTAATGGACCGAGCGAAGGAGAGCTCTCAGGGACAGATACTAACTCTCAGAGACTGGGACTA  
CTACTCAATCAGCGGGCCTTGACCACACTACTAAACCACCGGGAGGTCAAGTCTCCACGAAGCCACCAGC  
TGTGACCAGGGAAGGCTTAGGTGTCGAGAAAAGCAGGGCACCTGCCCCAGCGTGGACATACCCAAGCTC  
GGCCTCTGTGAGGACCAGTGTGAGGTGGACAGCCAGTGTCTGGCAACATGAAATGCTGCCGCAATGGAT  
GTGGGAAGATGGCCTGCACCACACCCAAATTC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR201500 protein sequence  
Red=Cloning site Green=Tags(s)

MPACRLCLLAAGLLLGLLLFTPISATGTDAEKPGCEPQLEPITDCVLECTLDKDCADNRKCCQAGCSSVC  
SKPNGPSEGELSGTDTKLSETGTTTQAGLDHTTKPPGGQVSTKPPAVTREGLVREKQGTCPSPVDIPKL  
GLCEDQCQVDSQCSGNMKCCRNCGKMACTTPKF

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** BC099427

**ORF Size:** 522 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC099427](#), [AAH99427](#)

**RefSeq Size:** 737 bp

**RefSeq ORF:** 524 bp


**Locus ID:** 67701

**Cytogenetics:** 2 H3

**MW:** 18 kDa

Gene Summary: Broad range protease inhibitor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201500