

Product datasheet for **MR201505L4V**

Cd99 (NM_025584) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd99 (NM_025584) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd99
Synonyms:	1110061M03Rik; 2410026K10Rik; D4; pilr-1; Pilr-l
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025584
ORF Size:	525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201505).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025584.2 , NP_079860.2
RefSeq Size:	950 bp
RefSeq ORF:	528 bp
Locus ID:	673094
UniProt ID:	Q8VCN6
Cytogenetics:	4
Gene Summary:	Involved in T-cell adhesion processes. Plays a role in a late step of leukocyte extravasation helping leukocytes to overcome the endothelial basement membrane. Acts at the same site as, but independently of, PECAM1.[UniProtKB/Swiss-Prot Function]



[View online »](#)