

Product datasheet for **MR201506**

Myl7 (NM_022879) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Myl7 (NM_022879) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Myl7
Synonyms: MLC-2alpha; MLC2a; MYL2A; Mylc2a; RLC-A
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201506 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAGTAGGAAGGCTGGGACCCGAGGCAAGGCTGCAGCCACCAAGCAGGCACAACGTGGCTCTTCTA
 ATGTCTTCTCAATGTTCGAACAAGCCCAGATTCAGGAATCAAGGAAGCCTTCAGCTGCATTGACCAGAA
 CAGGGATGGGATCATCTGCAATCAGACCTGAAGGAGACCTATCCAGCTCGGGAGGGTAAGTGTCCG
 GAGGAAGAGCTGGACCCATGCTGCAGGAAGGGAAGGGTCCCATCAACTCACCGTCTTCTCACACTCT
 TCGGGGAGAAGCTCAATGGGACGACCCCGAGGAAGCCATCCTGAGTGCCTTCCGCATGTTGACCCAG
 CGGCCAGGGGTGGTGAACAAGGAAGAGTTCAAGCAGCTTCTCATGACCCAGGCAGACAAGTTCTCTCT
 GCTGAGGTAGAGCAACTGTTTGCCTGACACCCATGGACCTGGCCGCAACATTGACTACAAGTCTCTCT
 GCTACATCATACCCACGGGGATGAGAAAGAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201506 protein sequence
 Red=Cloning site Green=Tags(s)

MASRKAGTRGKAAATKQAQRGSSNVFSMFEQAQIQEFKEAFSCIDQNRDGIICKSDLKETYSQLGRVSV
 EEELDAMLQEGKGPINFTVFLTLFGEKLNKTDPEEAILSAFRMFDPSPGQGVVNKEEFKQLLMTQADKFS
 AEVEQLFALTPMDLAGNIDYKSLCYIITHGDEKEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_022879

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022879.2](#), [NP_075017.2](#)

RefSeq Size: 599 bp

RefSeq ORF: 528 bp


Locus ID: 17898

UniProt ID: [Q9QVP4](#)

Cytogenetics: 11 3.88 cM

MW: 19.4 kDa

Product images:



Circular map for MR201506