

Product datasheet for **MR201509L3V**

Wrb (NM_207301) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wrb (NM_207301) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wrb
Synonyms:	5530402J05Rik; C030018G21Rik; Chd5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207301
ORF Size:	525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201509).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207301.2
RefSeq Size:	2483 bp
RefSeq ORF:	525 bp
Locus ID:	71446
UniProt ID:	Q8K0D7
Cytogenetics:	16 C4
Gene Summary:	Receptor for ASNA1/TRC40-mediated insertion of tail-anchored (TA) proteins into the ER membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)