

Product datasheet for MR201523

Idnk (NM_198004) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Idnk (NM_198004) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Idnk
Synonyms:	5133401N09Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR201523 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGTGAGCGGCTCCGAAAATCCACTGTGGGTGCGCTGCTGGCCTCCAAGCTGGGATGGAAATTC
ATGATGCCGATGATTACCACTCCGAGGAGAATCGGATAAAGATGGCGAAAGGGTACCCTGAGCGACCA
GGACAGGATCCATGGCTTTCACCTTGCACGACATTTTACTAAGAGATGTGGCTTTGGGACAGCCTGTT
GTTCTAGCCTGTTAGCTCTGAAGAAAACGTACAGAGACATCTTGATCCGAGGAGGAAGTATGCGCCCC
TGAAAAGTGATGACTCAGCAAAGGAACCGCTGGCTGGTGGGAAGCTCCTGGTGGTCTACCTGCGGCTC
GTTTGACATCATTATGGACGCTTGCTCCAAGAAAAGGACATTTTATGCCGCCGAGTTACTGCAGTCC
CAGTTTAGTATTCTGGAGCCCCATCAGCTCCCGAAAACCTTCATCCAAGTCAGCGTGACAAAAGTCTCC
CGGAGATCACCGCTGCCGTTATGGAAGCCCTCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR201523 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MGVSGSGKSTVGALLASKLGWKFYDADDYHSEENRIKMAKGVPLSDQDRIPWLCTLDILLRDVALGQPV
VLACSALKKTYRDILIRGGSDAPLKSDDSAKEPLAGKLLVVYLCGSFDIYGRLLQRKGHFMPPELLQS
QFSILEPPSAPENFIQVSVDKSLPEITAAVMEALK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	Sgfl-MluI
--------------------	-----------



[View online »](#)

Cloning Scheme:


ACCN: NM_198004

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_198004.1](#)

RefSeq Size: 1511 bp

RefSeq ORF: 555 bp


Locus ID: 75731

UniProt ID: [Q8R0J8](#)

Cytogenetics: 13 B1

MW: 19.1 kDa

Product images:



Circular map for MR201523