

## Product datasheet for **MR201535L3V**

### Tmem97 (NM\_133706) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem97 (NM_133706) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem97
Synonyms:	1810014L12Rik; AI115531; AL022956; D11Bhm182e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133706
ORF Size:	531 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133706.2</a> , <a href="#">NP_598467.1</a>
RefSeq Size:	1333 bp
RefSeq ORF:	531 bp
Locus ID:	69071
UniProt ID:	<a href="#">Q8VD00</a>
Cytogenetics:	11 46.74 cM



[View online »](#)

**Gene Summary:**

Intracellular orphan receptor that binds numerous drugs and which is highly expressed in various proliferating cells. Corresponds to the sigma-2 receptor, which is thought to play important role in regulating cell survival, morphology and differentiation. May play a role as a regulator of cellular cholesterol homeostasis. May function as sterol isomerase. May alter the activity of some cytochrome P450 proteins.[UniProtKB/Swiss-Prot Function]