

Product datasheet for MR201544

Ndufb10 (NM_026684) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ndufb10 (NM_026684) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ndufb10
Synonyms: 22kDa; 0610011B04Rik; PDSW
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201544 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCAGATTCTGGGACAAGGATGTGTACCCGGAGCCCCGAGCCGCACGCCTGCTCCCTCGCCGAGA
 CCTCGCTCCCTAACCCCATCACCTACTTGACGAAGGCCTACGACCTCGTCGTGGACTGGCCCGTGACCT
 CGTGAGAGAGTTTATTGAACGACAGCATGCCAAGAACCGAACCTACTACTACCACCGACAGTACCGTCGA
 GTGCCAGACATCAGAAATGCAAAGAGGGTGATGTTCTGTGTATCTATGAGGCTGAGATGCAGTGGAGAA
 GGGACTCAAAGTTGACCAGGAGATCATGAACATCATCCAGGAGAGACTGAAGGCTTGTCAGCAGAGGGA
 GGGAGAGAACTACCAGCAGAACTGTGCCAAGGAACTGGAGCAGTTCACCAAAGTGACTAAAGCCTACCAG
 GACCGCTACTTGACCTGGGAGCCTATTATTCTGCCAGAAAGTGTCTGGCAAAGCAGAAGCAGAGAATGC
 TGAAGAAAGGAAGGCTGCGAGGAGGAGGCTGCTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201544 protein sequence
 Red=Cloning site Green=Tags(s)

MPDSWDKDVYPEPPSRTPAPSPQTSPLPNPITYLTKAYDLVVDWPVTLVREFIERQHAKNRTYYYHRQYRR
 VPDITECKEGDVLCIYEAEQWRRDFKVDQEIMNIIQERLKACQQREGENYQQNCAKELEQFTKVTKAYQ
 DRYLDLGAYYSARKCLAKQKQRMLEERKAARQEAAA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:



ACCN: NM_026684

ORF Size: 531 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026684.2](#)

RefSeq Size: 657 bp

RefSeq ORF: 531 bp

Locus ID: 68342


UniProt ID: [Q9DCS9](#)

Cytogenetics: 17 A3.3

MW: 21 kDa

Gene Summary: Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201544