

Product datasheet for **MR201555L4V**

Odf4 (BC060992) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Odf4 (BC060992) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Odf4 |
| Synonyms: | Oppo1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | BC060992 |
| ORF Size: | 531 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR201555). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | BC060992.1 |
| RefSeq Size: | 1332 bp |
| RefSeq ORF: | 533 bp |
| Locus ID: | 252868 |
| Cytogenetics: | 11 B3 |
| Gene Summary: | Component of the outer dense fibers (ODF) of spermatozoa which could be involved in sperm tail structure, sperm movement and general organization of cellular cytoskeleton. [UniProtKB/Swiss-Prot Function] |



[View online »](#)