

Product datasheet for **MR201571**

3110082I17Rik (BC058580) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 3110082I17Rik (BC058580) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 3110082I17Rik
Synonyms: AW411904
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201571 representing BC058580
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAACATAAGCGAAAAGGTCTCGAAGGGACAGGAAAGGAGAGTAAAAGGCAGAAAATAACACCGG
CAGAGGAGACGCCGAGGACATCGGAGGCTGGGCCAGACAAGGAGACGCCAGTACCCTGGTGCAGGAGGC
CAGCCCTGAGTTGTCCCCAGAGGAGAGGAGAGTCTTGAAAGGAACTGAAGAAGGAACGGAAGAAGGAA
GAGAAGAAACGGCTGAGGAGGGCGGCATCGCTGCTACACAGACTGCCAAGGTTCCAGACGCTGCCAGCCA
AGCCCTCGGCTGCCACCCTGGCTCTGGAGTACCTCCAAGGGTGGGCTCAGAAGCAGGAGAGCTGGCGGTT
CCAGAAGACAAGACAGACGTGGCTCCTGCTGCACATGTATGATGAGGACAAGGTTCTGTGCACATGGATGG
GCCAAGGGGGGCTGGGGGGCCGTGGGACATTCAGGTCACAGGGTCTCGGCATGGCAGGTGGGCACCTTGGG
TGTCCGGGTCCATGCACAGCACATGGACAACCTAAGCAAAGTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201571 representing BC058580
Red=Cloning site Green=Tags(s)
MAKHKRKGLEGTGKESKRQKITPAEETPRTPSEAGPKETASTLVQEASPELSPEERRVLERKLLKKERKKE
EKKRLREAGIAATQTAKVQTLPAKPSAATLALAYLQGWAKQESWRFQKTRQTWLLLHMYDEDKVCAGHW
AKGGWGA VGHSGHRVSAWQVGTWVSGSMHSTWTTKQSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: BC058580

ORF Size: 534 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC058580](#)

RefSeq Size: 2059 bp


RefSeq ORF: 536 bp

Locus ID: 73212

Cytogenetics: 5 G2

MW: 75.5 kDa

Product images:



Circular map for MR201571