

Product datasheet for **MR201592**

Arl5c (NM_207231) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arl5c (NM_207231) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Arl5c
Synonyms: Arl12; BC065791
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201592 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGACAGCTGATAGCCAAGCTGATGCGAATCTTTGGGAGTCAGGAGCACAAAGTCATTATCGTGGGAC
 TGGACAATGCAGGAAGACCACCATCTCTACCAGTTCCTGACCAATGAGGTTGTTCACACCTGCTCCAC
 CATTGGGAGCAACGTGGAAGAAATTGTTCTTCGGAAGACTCACTTCCTCATGTGGGACCTAGGGGGCCAG
 GAGGCGCTGCGCTCCACCTGGGACACTTACTACTCTAATGCTGAGTTTGTATCCTTGTGATAGATAGCA
 CAGACCGGAATCGGCTGCTCACCACCGGGAGGAGCTGTATAAGATGCTGGCCCATGAGGCACCTCAAAA
 TGCTTCAGTACTGATCTTTGCCAACAAGCAGGATGTGAAAGACTCCATGACCACAGCAGAGATCTCCAA
 TTTCTCACTCAGTGCCATCAAAGACCACCGTGGCATATCCAGGGCTGCTGTGCTCTCACAGGAGAAG
 GGCTGCCAGCTGGGCTCCAGTGGATGCAAGCTCAGGCCACCGCCAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201592 protein sequence
 Red=Cloning site Green=Tags(s)


MGQLIAKLMRIFGSQEHKVIIVGLDNAGKTTILYQFLTNEVVHTCSTIGSNVEEIVLRKTHFLMWDLGGQ
 EALRSTWDITYSNAEFVILVIDSTDRNRLTTREELYKMLAHEALQNASVLI FANKQDVKDSMTTAEISQ
 FLTLSAIKDHPWHIQGCCALTGEGLPAGLQWMQAQATAN

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)


Cloning Scheme:



ACCN:	NM_207231
ORF Size:	537 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_207231.1, NP_997114.1</u>
RefSeq Size:	1616 bp
RefSeq ORF:	540 bp
Locus ID:	217151

UniProt ID: [Q6P068](#)
 Cytogenetics: 11 D
 MW: 20.1 kDa
 Gene Summary: Binds and exchanges GTP and GDP.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201592