

## Product datasheet for MR201598

### Yaf2 (NM\_024189) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Yaf2 (NM\_024189) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Yaf2  
**Synonyms:** 2810021M11Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201598 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGCGACAAGAAGAGCCCCACCAGGCCGAAGCGGCAGCCGAAGCCCGCCTCGGATGAGGGCTACTGGG  
 ACTGCAGTGTCTGCACCTTCCGGAACAGCGCCGAGGCCTCAAGTGCATGATGTGCGATGTGCGGAAGGG  
 CACCTCCACCCGGAAGCCACGCCCTGTCTCCAGCTGGTGCACAGCAGGTCACGCAGCAGTTCTGTGCC  
 CCTACCCAGTCGAAGAAAGAGAAGAAAGACCGAGTAGAGAAGGACAAGAGTGAGAAGGAAGCTGCTAGCA  
 AAAAGAAGTGTACAAGAAAACCAGACCAAGATTAATAAATGTGGATCGGAGTAGCGCTCAGCATTGGA  
 AGTCACTGTTGGAGATCTAACAGTCATTATCACAGACTTTAAGGAGAAAGCCAAGTCAGCGCCTGCGTCC  
 AGCGCTGCTGGCGATCAGCACAGCCAGGGCAGCTGCAGCTCGGACAGCACTGAGCGGGCGTGTCCAGGT  
 CCTCTCGCCGCGGGGAGAAGCCTCCTCCCTCAACGGAGAGTCTCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR201598 protein sequence  
 Red=Cloning site Green=Tags(s)

MGDKKSPTRPKRQPKPASDEGYWDCSVCTFRNSAEAFKCMCDVRKGTSTRKPRPVSQQLVAQQVTQQFVP  
 PTQSKKEKDRVEKDKSEKEAASKKNCHKTRPRLKNVDRSSAQHLEVTVGDLTVIITDFKEKAKSAPAS  
 SAAGDQHSQGSSSDSTERGVSRSSSPRGEASSLNGESH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NM\_024189

**ORF Size:** 540 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024189.1](#)

**RefSeq Size:** 2248 bp

**RefSeq ORF:** 540 bp

**Locus ID:** 67057


**UniProt ID:** [Q99LW6](#)

**Cytogenetics:** 15 E3

**MW:** 19.7 kDa

**Gene Summary:** Binds to MYC and inhibits MYC-mediated transactivation. Also binds to MYCN and enhances MYCN-dependent transcriptional activation. Increases calpain 2-mediated proteolysis of YY1 in vitro. Component of the E2F6.com-1 complex, a repressive complex that methylates 'Lys-9' of histone H3, suggesting that it is involved in chromatin-remodeling.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR201598