

Product datasheet for **MR201615A1V**

Mouse Ccdc115 (NM_027159) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccdc115 (NM_027159) AAV Particle
Tag: Myc-DDK
Symbol: Ccdc115
Synonyms: 2310061I09Rik; Ccp1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201615 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGTCCAGGCTCTGCGAGAGGAGTTGGACTCCAAGTGCCTACAGCTGCTCAGTGACCTGGAGGAAC
 TGGAGGCGAAGCGGGCGCGCTGAACGCCAGGGTGGAGGAGGGTTGGCTCTCACTTGCCAAGGCTCGCTA
 TGCCATGGGAGCCAAGTCGGTAGGACCCCTGCAGTATGCCTCGCGTATGGAGCCTCAGGTCTGCGTGCCT
 GCCAGCGAGGCCCAGGATGGACCCAGACCTTCAGGGTGATCAAAGCTGACGCCAGACCCCGAGGAAG
 TGGGGCCTAGCGAAGCATCTCTGCGCAGGCGCAAGGGCCCAAGACCAAGAGTTAGGATCTGCTGT
 AGTTCCCTCAAGATCCCTTGAAGTGGTTTGGAAATTTGGTTCTCACAGTCTGCGGCAGGCCCAAGCCAGC
 TTCCGGGATGGCTACAGCTGGCTGCAGATATAGCCAGTCTCCAGACTCGCATCAACTGGGGTCAGAGTC
 AGCTTCGAGGCCTCCAGAAAAAAGTGAAGGAGTTGGACCCTGGGCCTGCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201615 protein sequence
 Red=Cloning site Green=Tags(s)

MAVQALREELDSKCLQLLSDLEELEAKRAALNARVEEGWLSLAKARYAMGAKSVGPLQYASRMEPQVCVR
 ASEAQDGPQTRVIKADAQTPEEVGPSEASLRRRKGPTKTKELGSAVVPQDPLNWFILVPHSLRQAQAS
 FRDGLQLAADIASLQTRINWQSQLRGLQKKLKELDPGPA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2



ACCN:	NM_027159
ORF Size:	540 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_027159.1</u>
RefSeq Size:	1721 bp
RefSeq ORF:	543 bp
Locus ID:	69668
UniProt ID:	<u>Q8VE99</u>
Cytogenetics:	1 B
MW:	19.7 kDa