

Product datasheet for **MR201619L3V**

II17b (NM_019508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	II17b (NM_019508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	II17b
Synonyms:	1110006O16Rik; 1700006N07Rik; Zcyto7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019508
ORF Size:	543 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201619).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019508.1 , NP_062381.1
RefSeq Size:	692 bp
RefSeq ORF:	543 bp
Locus ID:	56069
UniProt ID:	Q9QXT6
Cytogenetics:	18 E1
Gene Summary:	Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.[UniProtKB/Swiss-Prot Function]



[View online »](#)