

Product datasheet for **MR201629L4V**

Arf3 (NM_007478) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arf3 (NM_007478) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arf3
Synonyms:	5430400P17Rik; AI854770
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007478
ORF Size:	546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201629).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007478.3 , NP_031504.1
RefSeq Size:	3481 bp
RefSeq ORF:	546 bp
Locus ID:	11842
UniProt ID:	P61205
Cytogenetics:	15 F1
Gene Summary:	GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.[UniProtKB/Swiss-Prot Function]



[View online »](#)