

Product datasheet for **MR201633A1V**

Mouse Csn3 (NM_007786) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Csn3 (NM_007786) AAV Particle
Tag: Myc-DDK
Symbol: Csn3
Synonyms: AW208918; CSN10; Csnk
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201633 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGAGGAATTTTATCGTAGTTGTGAATATTCTGGCATTAACTCTGCCCTTTTGGCTGCAGAGATAC
 AAAACCCAGATTCAAACCTGCCGTGGTGAGAAGAATGACATAGTGTACGATGAGCAAAGAGTCCTGTACAC
 ACCAGTCCGCTCCGTCTGAACCTCAATCAATATGAACCTAATTACTACCACTACAGGCCATCTTTGCCCT
 GCTACTGCTAGCCCATACATGTACTACCACTTGTGTAGACTACTTCTGCTTAGGTCTCCAGCCCAA
 TCTCCAAATGGCAATCAATGCCAAATTTCCCCAGTCTGCTGGAGTACCTTATGCCATTCCAAACCCATC
 CTTTCTTGCCATGCCTACAAATGAAATCAGGATAACACCGCCATTCCAACCATCGACCCATTACTCCC
 ATTGTGTCCACACCAGTTCCTACCATGGAGTCAATAGTAAATACTGTGGCCAATCCTGAAGCTTCCACAG
 TGTCTATCAATACCCCTGAGACTACTACAGTCCCACTGCTCAACTGCAGCA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201633 protein sequence
 Red=Cloning site Green=Tags(s)

MMRNFIVVNI LALTLPLAAEIQNPD SNCRGEKNDIVYDEQ RVLVTPVRSVLNFNQYEPNYYHYRPSLP
 ATASPYMYPLVVRLLLLRSPAPISKWQSMNFPQSAGVPYAIPNPSFLAMPTNENQDNTAIPITP
 IVSTPVPTMESIVNTVANPEASTVSINTPETTTVPVSSTAA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2


[View online »](#)

ACCN:	NM_007786
ORF Size:	543 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_007786.2
RefSeq Size:	1597 bp
RefSeq ORF:	546 bp
Locus ID:	12994
UniProt ID:	P06796
Cytogenetics:	5 43.56 cM
MW:	20 kDa