

Product datasheet for **MR201647A1V**

Mouse Mid1ip1 (NM_026524) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Mid1ip1 (NM_026524) AAV Particle
Tag: Myc-DDK
Symbol: Mid1ip1
Synonyms: 3110038L01Rik; Mig12; S14R
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201647 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGCAAACTCTGCGACACATATAACCAGAAGCACTCGCTCTTTAACGCCATGAATCGCTTCATTGGCG
 CGGTGAACAACATGGACCAGACGGTGATGGTGCCAGTCTGCTGCGCGACGTACCCCTGTCCGAGCCTGA
 GATAGACGAGGTGAGCGTGGAGGTAGGCGGCAGTGGCGGCTGCCTGGAGGAGCGCAGACCCCGGCCCA
 AGCCCGGGCAGCGCAATGAAAGCTTTTCGCGCCCTCCCGGGACATGTACAGCCACTACGTGCTCCTCA
 AGTCCATCCGCAATGATATCGAGTGGGAGTCTGCACCAGCCTTCGTCTCCGCCGGCCGGGAGCGAGGA
 GAGCACCTGGAAGCCCAAGGACATCCTGGTGGCCTGAGTCACTTGGAGAGCGCGGATGCGGGCGAGGAA
 GATCTGGAGCAGCAGTTCCACTACCACCTGCGCGGGCTGCACACCGTGTCTCCAAACTACCCGAAAAG
 CCAACATCCTCACCAATAGATACAAGCAGGAGATCGGCTTCAGTAATTGGGGCCAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201647 protein sequence
 Red=Cloning site Green=Tags(s)

MMQICDTYNQKHSLFNAMNRFIGAVNNMDQTMVPSLLRDVPLSEPEIDEVSVEVGGSGGCLEERTTPAP
 SPGSANESFFAPSRDMYSHYVLLKSIRNDIEWGLVHQPSSPPAGSEESTWPKDILVGLSHLESADAGEE
 DLEQQFHYHLRGLHTVLSKLTRKANILTNRKQEIGFSNWH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2


[View online »](#)

ACCN:	NM_026524
ORF Size:	546 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026524.2</u>
RefSeq Size:	2056 bp
RefSeq ORF:	549 bp
Locus ID:	68041
UniProt ID:	<u>Q9CQ20</u>
Cytogenetics:	X A1.1
MW:	20.4 kDa