

Product datasheet for **MR201647L4V**

Mid1ip1 (NM_026524) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mid1ip1 (NM_026524) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mid1ip1
Synonyms:	3110038L01Rik; Mig12; S14R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026524
ORF Size:	549 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201647).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026524.2
RefSeq Size:	2056 bp
RefSeq ORF:	549 bp
Locus ID:	68041
UniProt ID:	Q9CQ20
Cytogenetics:	X A1.1
Gene Summary:	Plays a role in the regulation of lipogenesis in liver. Up-regulates ACACA enzyme activity. Required for efficient lipid biosynthesis, including triacylglycerol, diacylglycerol and phospholipid. Involved in stabilization of microtubules.[UniProtKB/Swiss-Prot Function]



[View online »](#)