

Product datasheet for **MR201653L3V**

Arl3 (NM_019718) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arl3 (NM_019718) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arl3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019718
ORF Size:	549 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201653).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019718.2
RefSeq Size:	838 bp
RefSeq ORF:	549 bp
Locus ID:	56350
UniProt ID:	Q9WUL7
Cytogenetics:	19 C3



[View online »](#)

Gene Summary:

Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP) (PubMed:18376416). Required for normal cytokinesis and cilia signaling. Required for targeting proteins to the cilium, including myristoylated NPHP3 and prenylated INPP5E. Targets NPHP3 to the ciliary membrane by releasing myristoylated NPHP3 from UNC119B cargo adapter into the cilium (By similarity). Requires assistance from GTPase-activating proteins (GAPs) like RP2 and PDE6D, in order to cycle between inactive GDP-bound and active GTP-bound forms (PubMed:15979089). Required for PKD1:PKD2 complex targeting from the trans-Golgi network to the cilium (PubMed:25405894).[UniProtKB/Swiss-Prot Function]