

## Product datasheet for **MR201667L3V**

### Rap2b (NM\_028712) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rap2b (NM_028712) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rap2b
Synonyms:	4021402C18Rik; AA408554
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028712
ORF Size:	552 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028712.4</a> , <a href="#">NP_082988.1</a>
RefSeq Size:	4197 bp
RefSeq ORF:	552 bp
Locus ID:	74012
UniProt ID:	<a href="#">P61226</a>
Cytogenetics:	3 E1



[View online »](#)

**Gene Summary:**

Small GTP-binding protein which cycles between a GDP-bound inactive and a GTP-bound active form. Involved in EGFR and CHRM3 signaling pathways through stimulation of PLCE1. May play a role in cytoskeletal rearrangements and regulate cell spreading through activation of the effector TNIK. May regulate membrane vesiculation in red blood cells.  
[UniProtKB/Swiss-Prot Function]