

## Product datasheet for MR201688

### Mbp (NM\_001025254) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Mbp (NM_001025254) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mbp
Synonyms:	C76307; goll; golli-mbp; Hmb; Hmbpr; jv; jve; mld; R75289; shi
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR201688 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGGCATCACAGAAGAGACCCCTCACAGCGATCCAAGTACCTGGCCACAGCAAGTACCATGGACCATGCCA  
 GGCATGGCTTCTCCCAAGGCACAGAGACCGGCATCCTTGACTCCATCGGGCGCTTCTTAGCGGTGA  
 CAGGGGTGCGCCCAAGCGGGCTCTGGCAAGGTACCCTGGCTAAAGCAGAGCCGGAGCCCTCTGCCCTCT  
 CATGCCCGCAGCCGTCTGGGCTGTGCCACATGTACAAGGACTCACACACGAGAACTACCCATTATGGCT  
 CCCTGCCCCAGAAGTCGACGACGCGCCGACCCAAGATGAAAACCCAGTAGTCCATTTCTCAAGAACAT  
 TGTGACACCTCGAACACCACCTCCATCCCAAGGGAAGGGGCGGAGGGGCAGAAGCCAGGATTTGGCTAC  
 GGAGGCAGAGCTTCCGACTATAAATCGGCTCACAAGGGATTCAAGGGGGCCTACGACGCCAGGGACGC  
 TTTCCAAAATCTTTAAGCTGGGAGGAAGAGACGCGCTCTGGATCTCCCATGGCGAGACGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR201688 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MASQKRPSQRSKYLATASTMDHARHGFLPRHRDTGILDSIGRFFSGDRGAPKRGSGKVPWLKQSRSPGPS  
 HARSRPLCHMYKDSHTRTTHYGLPQKSQHGRTQDENPVVHFFKNIIVTPRTPPPSQKGAGQKPGFGY  
 GGRASDYKSAHKGFKGAYDAQGTLSKIFKLGGDRSRSGSPMARR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	Sgfl-MluI
--------------------	-----------


[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001025254

**ORF Size:** 552 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_001025254.2](#), [NP\\_001020425.1](#)

**RefSeq Size:** 2125 bp

**RefSeq ORF:** 555 bp

**Locus ID:** 17196


UniProt ID: [P04370](#)

Cytogenetics: 18 55.84 cM

MW: 20.3 kDa

**Gene Summary:** The protein encoded by the classic Mbp gene is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the nervous system. However, Mbp-related transcripts are also present in the bone marrow and the immune system. These mRNAs arise from the long Mbp gene (otherwise called "Golli-Mbp") that contains 3 additional exons located upstream of the classic Mbp exons. Alternative splicing from the Golli and the Mbp transcription start sites gives rise to 2 sets of Mbp-related transcripts and gene products. The Golli mRNAs contain 3 exons unique to Golli-Mbp, spliced in-frame to 1 or more Mbp exons. They encode hybrid proteins that have N-terminal Golli aa sequence linked to Mbp aa sequence. The second family of transcripts contain only Mbp exons and produce the well characterized myelin basic proteins. This complex gene structure is conserved among species suggesting that the Mbp transcription unit is an integral part of the Golli transcription unit and that this arrangement is important for the function and/or regulation of these genes. Mutation of the Mbp gene is associated with the 'shiverer' and 'myelin deficient' phenotypes in mouse. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for MR201688