

Product datasheet for **MR201693L4V**

Adm (NM_009627) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adm (NM_009627) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adm
Synonyms:	AM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009627
ORF Size:	555 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201693).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009627.1
RefSeq Size:	1381 bp
RefSeq ORF:	555 bp
Locus ID:	11535
UniProt ID:	P97297
Cytogenetics:	7 57.7 cM
Gene Summary:	AM and PAMP are potent hypotensive and vasodilator agents.[UniProtKB/Swiss-Prot Function]



[View online »](#)