

Product datasheet for MR201695

Gnpnat1 (NM_019425) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gnpnat1 (NM_019425) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gnpnat1
Synonyms: AU017428; AU040593; EMeg32; Gpnat1; Gsnpat
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201695 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAACCTGATGAAACTCCCATGTTTGACCCAAGTTTGCTCAAAGAAGTGGACTGGAGTCAGAATACAG
 CTATTTTTTCTCCAGCCATTTCTCCAACGCATCCTGGAGAAGGCTTGGTTTTGAGCCTCTCTGTACTGC
 AGACTTAAATAAAGGTTTTTTAAAGGATTGGGTACGCTGACAGAGACTGGCGTCGTCAGCCCGGAACAG
 TTCATGAAATCTTTGAGCATATGAAGAAGTCTGGAGATTACTACGTTACAGTTGTGGAAGATGTGACCC
 TAGGACAAATTGGGCTACAGCAACTCTGATAATAGAACATAAAATTTATCCATTCATGTGCTAAGAGAGG
 GAGAGTAGAAGATGTCGTCGTTAGTGACGAGTGCAGAGGGAAGCAACTCGGCAAACCTGTTATTATCAACT
 CTTACTTTGCTAAGCAAGAAATGAACTGTTACAAGATCACCCCTTGAATGTCTACCACAAAATGTTGGTT
 TCTACAAAAGTTTGACTATACAGTATCTGAAGAAAATTACATGTGTCGGAGGTTCTAAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201695 protein sequence
 Red=Cloning site Green=Tags(s)

MKPDETPMFDPSLLKEVDWSQNTAIFSPAISPTHPEGLVLRPLCTADLNKGFFKVLGQLTETGVVSPEQ
 FMKSFHEMKKSGDYVTVVVDVTLGQIVATATLIEHKFIHSCAKRGRVEDVVVSDECRGKQLGKLLLSL
 LTLLSKKLNKYKITLECLPQNVGFYKFDYTVSEENYMCRRFLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_019425

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019425.1](#), [NM_019425.2](#), [NP_062298.1](#)

RefSeq Size: 2536 bp

RefSeq ORF: 555 bp


Locus ID: 54342

UniProt ID: [Q9JK38](#)

Cytogenetics: 14 C1

MW: 20.8 kDa

Product images:



Circular map for MR201695