

Product datasheet for **MR201698L4V**

Tnfaip8l2 (NM_027206) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Tnfaip8l2 (NM_027206) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Tnfaip8l2 |
| Synonyms: | 1810019A08Rik; AW610835; TIPE2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_027206 |
| ORF Size: | 555 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR201698). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_027206.1 |
| RefSeq Size: | 1122 bp |
| RefSeq ORF: | 555 bp |
| Locus ID: | 69769 |
| UniProt ID: | Q9D8Y7 |
| Cytogenetics: | 3 F2.1 |



[View online »](#)

Gene Summary:

Acts as a negative regulator of innate and adaptive immunity by maintaining immune homeostasis. Negative regulator of Toll-like receptor and T-cell receptor function. Prevents hyperresponsiveness of the immune system and maintains immune homeostasis. Inhibits JUN/AP1 and NF-kappa-B activation. Promotes Fas-induced apoptosis.[UniProtKB/Swiss-Prot Function]