

Product datasheet for MR201716

Tnfrsf17 (NM_011608) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tnfrsf17 (NM_011608) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tnfrsf17
Synonyms: BCM; BCMA; Tnfrsf13; Tnfrsf13a
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201716 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGCAACAGTGTTCACAGTGAATATTTTGACAGTCTGCTGCATGCTTGCAAACCGTGCACTTGC
 GATGTTCCAACCTCCTGCAACCTGTCAGCCTTACTGTGATCCAAGCGTGACCAGTTCAGTGAAAGGGAC
 GTACACGGTGTCTGGATCTTCTGGGGCTGACCTTGGTCCTCTCTTTGGCACTTTTACAATCTCATT
 TTGCTGAGGAAGATGAACCCGAGGCCCTGAAGGACGAGCCTCAAAGCCAGGTCAGCTTGACGGATCGG
 CTCAGCTGGACAAGGCCGACACCGAGCTGACTAGGATCAGGGCTGGTGACGACAGGATCTTCCCGAAG
 CCTGGAGTATACAGTGAAGAGTGCACCTGTGAGGACTGTGTCAAGAGCAAACCAAGGGGATTCTGAC
 CATTCTTCCCGCTTCCAGCCATGGAGGAGGGGCAACCATTTGTCAACCACAAAAACGGGTGACTACG
 GCAAGTCAAGTGTCCAAGTCTTTGCAAAGTGTGATGGGGATGGAGAAGCCAACCTCACACTAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201716 protein sequence
 Red=Cloning site Green=Tags(s)

MAQQCFHSEYFDSLLHACKPCHLRCSNPPATCQPYCDPSVTSSVKGYTVLWIFLGLTLVLSLALFTISF
 LLRKMNPEALKDEPQSPGQLDGSALDKADTEL TRIRAGDDRIFPRSLEYTVEECTCEDCVKSKPKGDS
 HFFPLPAMEEGATILVTTKTGDYKSSVPTALQSVMGMEKPTHTR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_011608

ORF Size: 558 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011608.1](#), [NP_035738.1](#)

RefSeq Size:	812 bp
RefSeq ORF:	558 bp
Locus ID:	21935
UniProt ID:	<u>O88472</u>
Cytogenetics:	16 A1
MW:	20.4 kDa
Gene Summary:	Receptor for TNFSF13B/BLYS/BAFF and TNFSF13/APRIL. Promotes B-cell survival and plays a role in the regulation of humoral immunity. Activates NF-kappa-B and JNK (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201716