

## Product datasheet for **MR201716L1V**

### **Tnfrsf17 (NM\_011608) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Tnfrsf17 (NM_011608) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfrsf17
Synonyms:	BCM; BCMA; Tnfrsf13; Tnfrsf13a
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_011608
ORF Size:	558 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201716).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011608.1</a> , <a href="#">NP_035738.1</a>
RefSeq Size:	812 bp
RefSeq ORF:	558 bp
Locus ID:	21935
UniProt ID:	<a href="#">O88472</a>
Cytogenetics:	16 A1
Gene Summary:	Receptor for TNFSF13B/BLyS/BAFF and TNFSF13/APRIL. Promotes B-cell survival and plays a role in the regulation of humoral immunity. Activates NF-kappa-B and JNK (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)