

Product datasheet for **MR201717A1V**

Mouse Bmf (BC079650) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Bmf (BC079650) AAV Particle
Tag: Myc-DDK
Symbol: Bmf
Synonyms: AW260063; AW411872
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201717 representing BC079650
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCCACCTCAGTGTGTGGAGGAGCTAGAAGATGATGTGTTCCAGTCAGAGGATGGGGAGCCAGGGA
CACAGCCTGGGGCTTGCTCTCTGCTGACCTGTTTGCCAGAGCCAGCTGGACTGTCCCCTCAGTCGACT
CCAGCTCTTCCCTCTCACCCATTGCTGTGGTCCCGGACTCCGGCCCATAGCCAGGAAGACAAGGCCACT
CAGACCCTCAGTCCAGCTTCCCCAAGCCAGGGTGTCTGCTGCCTTGTGGGTGACAGAGGAACCCAGAG
GACTCTTTTACGGCAACGCTGGCTACAGGCTTCCTCTCCCTGCCAGTTTCCCTGCAGGCTCACCCCTTGG
GGAGCAGCCCCCTGAAGGACAGTTCCTTCAGCACCAGCAGAGGTGCAGATCGCCAGAAAGCTTCAGTGT
ATTGCAGACCAGTTCATCGGCTTCATACGCAACAACACCAGCAGAACCGAGACCGTGCCTGGTGGCAGG
TCTTCCTTTTCTTCAAACCTCGCCCTGAACAGACAAGAAAACAGGGAAGGGGTGGGGCCCTGG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201717 representing BC079650
Red=Cloning site Green=Tags(s)

MEPPQCVEELEDVDFQSEDGEPGTQPGGLLSADLFAQSQLDCPLSRLQLFPLTHCCGPGLRPISQEDKAT
QTLSPASPSQGVMLPCGVTEEPQRLFYGNAGYRLPLPASFPAGSPLGEQPPGQFLQHRAEVQIARKLQC
IADQFHRLHTQQHQNRDRAWVQVFLFLQNLALNRQENREGVGPW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2



[View online »](#)

ACCN:	BC079650
ORF Size:	555 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC079650.1
RefSeq Size:	4667 bp
RefSeq ORF:	557 bp
Locus ID:	171543
Cytogenetics:	2 E5
MW:	171.1 kDa