

## Product datasheet for **MR201739L3V**

### Ift27 (NM\_025931) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ift27 (NM_025931) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ift27
Synonyms:	2600013G09Rik; Rabl4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025931
ORF Size:	561 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025931.2</a>
RefSeq Size:	1050 bp
RefSeq ORF:	561 bp
Locus ID:	67042
UniProt ID:	<a href="#">Q9D0P8</a>
Cytogenetics:	15 E1



[View online »](#)

**Gene Summary:**

Small GTPase-like component of the intraflagellar transport (IFT) complex B that promotes the exit of the BBSome complex from cilia via its interaction with ARL6 (PubMed:25446516). Not involved in entry of the BBSome complex into cilium. Prevents aggregation of GTP-free ARL6. Required for hedgehog signaling (PubMed:25446516). Forms a subcomplex within the IFT complex B with IFT25 (By similarity). Its role in intraflagellar transport is mainly seen in tissues rich in ciliated cells such as kidney and testis. Essential for male fertility, spermiogenesis and sperm flagella formation (PubMed:28964737). Plays a role in the early development of the kidney (PubMed:29626631). May be involved in the regulation of ureteric bud initiation (PubMed:29626631).[UniProtKB/Swiss-Prot Function]