

Product datasheet for MR201747

Dctn3 (NM_016890) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dctn3 (NM_016890) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dctn3
Synonyms: p24
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201747 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCTCTGACCGATGTGCAGCGGCTACAGTCTCGAGTGGAGGAGCTGGAGCGCTGGGTGTACGGC
CGGGCGGGACCCGAGGCTCCAGGAAGGTGGCTGACGGTCTGGTCAAGGTGCAGGTGGCTTAGGGAACAT
TGCCAGTAAGAGGGAGAGGGTGAAGATTCTCTACAAAAAGATTGAAGACCTGATTAATACTGGATCCT
GAGTACATTGATCGATTGCCATACCTGAGGCCTCGAAGTTCAGTTCATCTTAGCAGAGGAGCAGTTCA
TCCTCTCCAGGTGGCACTCCTGGAACAAGTGAACGCCCTGGTGCCTGTGCTGGACAGCGCTTCTATCAA
AGCTGTTCTGAGCATGCTGCCCGCTGCAGCGCCTGGCGCAGATCCACATCCAGCAGCAGGACAGTGT
GTGGCGATCACTGAAGAGTCTAAGGCTCTCTGGAAGGGTATAACAAGACTACAATGCTTCTCTCCAAGC
AGTTTGTGCAGTGGGATGAGTACTCTGCCAGCTAGAGGCTGCCAAGCAAGTGAAGCCAGCAGAGGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201747 protein sequence
Red=Cloning site Green=Tags(s)

MAALTDVQRLQSRVEELERWVYGPGGTRGSRKVADGLVKVQVALGNIASKRERVKILYKKIEDLIKYLDP
EYIDRIAIPASKLQFILAEEQFILSQVALLEQVNALVPVLDASAIKAVPEHAARLQRLAQIHIQQDQC
VAITEESKALLEGYNKTTMLLSKQFVQWDELLCQLEAAKQVKPAEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_016890

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016890.4](#), [NP_058586.3](#)

RefSeq Size: 966 bp

RefSeq ORF: 561 bp

Locus ID: 53598


UniProt ID: [Q9Z0Y1](#)

Cytogenetics: 4 A5

MW: 21 kDa

Gene Summary: Together with dynein may be involved in spindle assembly and cytokinesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201747