

Product datasheet for MR201753

Bcl2l12 (BC039807) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bcl2l12 (BC039807) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bcl2l12
Synonyms: Bcl-L12, Bcl2-L12
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201753 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGCCGGCTCTGAAGAGTTGGGTCTCCGAGAGGACACACTGAAGTTCTGACGGCCTTCCTTAAGCGCG
 GTGAAGTCGCTGGGTCCCCTGTCCCACTCCACCCAGGAGTCTGCCAAGAGGAGACAACAGATTTCT
 GAGCCGCCTCGAAGATGCCTCCCTGTCCCTGGGGAGAGGAGCTCCCCACGGAGTCTCTAGGCC
 CACTTCTTCCCCTCGCCCTGCTATGGTTCGAGCCTGGCCAGCCTTCCGAGTCTATGCCCTGG
 TTGCCAGCGGCTGGAACAATTGGTGCAGGAGCAACTGAAGTCCCGCCAGTTCAGAGTCCAGGGACC
 CCCACCCACGGAGAAGGAAGCCCTGCTAAGGCGGCTGGTAGCCTTGTGGAGGAAGAGGCAGAAGTCATC
 AACCAAAAGGTGACCAGCCTCTACTCCACGCTATGGAGAGGACAGTTGTGGAATAGGGGGCAGACGCTGG
 CGCTGCTTTGTTTGTGCTGGGTGTGGTTCAGGGCTCCAGTAAACATCAGAGCATCTTACAAGAGCC
 T

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201753 protein sequence
 Red=Cloning site Green=Tags(s)

MAGSEELGLREDTLKVLTAFLKRGEVAGSPVPTPPRSPAQEETDFLSRLRRCLPCPLGRGAPPTSSRP
 HFLPLRPCYGSEPGPASSEFYALVAQRLEQLVQEQLKSPSPSEFQGPPTTEKEALLRRLVALLLEEEAEVI
 NQKVTSLSYTLWRGQLWNRGQTLALLCFDAGCGSRAPVNIRASSQEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC039807

ORF Size: 561 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC039807](#), [AAH39807](#)

RefSeq Size: 1620 bp


RefSeq ORF: 563 bp

Locus ID: 75736

Cytogenetics: 7 B3

MW: 20.5 kDa

Product images:



Circular map for MR201753