

Product datasheet for **MR201753A1V**

Mouse Bcl2l12 (BC039807) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Bcl2l12 (BC039807) AAV Particle
Tag: Myc-DDK
Symbol: Bcl2l12
Synonyms: Bcl-L12, Bcl2-L12
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201753 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCGGCTCTGAAGAGTTGGGTCTCCGAGAGGACACACTGAAGGTTCTGACGGCCTTCCTTAAGCGCG
GTGAAGTCGCTGGGTCCCTGTCCCACTCCACCCAGGAGTCTGCCCAAGAGGAGACAACAGATTTCT
GAGCCGCTCCGAAGATGCCTCCCTGTCCCTGGGGAGAGGAGCTCCCCACGGAGTCCTTAGGCC
CACTTCCTTCCCTCCGCCCTGCTATGGTCCGAGCCTGGCCAGCCTCTTCCGAGTTCTATGCCCTGG
TTGCCAGCGGCTGGAACAATTGGTGCAGGAGCAACTGAAGTCCCCGCCAGTTCAGAGTTCAGGGACC
CCCACCCACGGAGAAGGAAGCCCTGCTAAGGCGGCTGGTAGCCTTGTGGAGGAAGAGGCAGAGTCATC
AACCAAAAGGTGACCAGCCTCTACTCCACGCTATGGAGAGGACAGTTGTGGAATAGGGGCGAGACGCTGG
CGCTGCTTTGTTTGTGCTGGGTGTGGTTCCAGGGCTCCAGTAAACATCAGAGCATCTTACAAGAGCC
T

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201753 protein sequence
Red=Cloning site Green=Tags(s)

MAGSEELGLREDTLKVLTAFLKRGEVAGSPVPTPPRSPAQEETDFLSRLRRCLPCPLGRGAPPTESSRP
HFLPLRPCYGSEPGPASSEFYALVAQRLEQLVQEQLKSPPSSEFQGGPPTEKEALLRRLVALLEEEAEVI
NQKVTSLYSTLWRGQLWNRGQTLALLCFDAGCGSRAPVNIRASSQEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse



Serotype:	AAV-2
ACCN:	BC039807
ORF Size:	561 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC039807</u> , <u>AAH39807</u>
RefSeq Size:	1620 bp
RefSeq ORF:	563 bp
Locus ID:	75736
Cytogenetics:	7 B3
MW:	20.5 kDa