

Product datasheet for MR201770

Pigh (BC086806) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pigh (BC086806) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pigh
Synonyms: 2210416H01Rik; A930028P05Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201770 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGACGAGAAGAGCTTCTCCGATATCTGCGCGGCCGCTAGCCCTGCGGTGCCGCTACTACTCAC
CGTACTGCCGCGAGTTCGGCCTCAGCTCGGCTCGGCTGTCGCTGTGCTCGCTCACC GCGGTACCTGCGC
GGTGTGGCTGGCGCCTACGGCTCTTACCCTGTGCGAGAATAGCATGGTCTCTCTGCCACCATCTC
ATCACCATCTTAGGCCTGCTTGGTTACCTCCATTTTGTGAAGATTGACCAGGAGACTCTACTAATCATCG
ACTCATTGGCATCCAGATGACTTCATCTTATGCATCGGGCAAAGAAAGTACCACCTTCATAGAGATGGA
CAAGGTCAAGGACATTATCATTAAACGAGGCTATCTACATGCAGAAGGTCACTTACTACCTCTGCATTTTA
CTGAAGGAGCCGGGAAGCCACACGAGATCTCTCGAGTAGTTCCTGTCTTCCAGAGTGCCAAAGCCCGGT
TGGACTGCTTGATTGAAGTATACAGAAGCTGCCAGGAGTTCTGGCACACCAGAAAGCCACAGCAACAAG
CTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201770 protein sequence
 Red=Cloning site Green=Tags(s)

MEDEKSFSDICGRLALRCRYSPYCREFLSSARLSLCSLTAVTCAVWLAAYGLFTLCENSMVLSATIF
ITILGLLGYLHFVKIDQETLLIIDSLGIQMTSSYASGKESTTFIEMDKVKDIIINEAIYMQKVIYYLCIL
LKEPGKPHEISRVPVVFQSAKPRLDCLIEVYRSCQEVLAHQKATATSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC086806

ORF Size: 564 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC086806](#), [AAH86806](#)

RefSeq Size: 2327 bp

RefSeq ORF: 566 bp


Locus ID: 110417

Cytogenetics: 12 35.51 cM

MW: 21.1 kDa

Gene Summary: Part of the complex catalyzing the transfer of N-acetylglucosamine from UDP-N-acetylglucosamine to phosphatidylinositol, the first step of GPI biosynthesis.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201770