

Product datasheet for MR201796L3V

Ramp2 (NM_019444) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ramp2 (NM_019444) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ramp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019444
ORF Size:	570 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201796).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019444.2 , NP_062317.1
RefSeq Size:	872 bp
RefSeq ORF:	570 bp
Locus ID:	54409
UniProt ID:	Q9WUP0
Cytogenetics:	11 64.47 cM
Gene Summary:	Transports the calcitonin gene-related peptide type 1 receptor (CALCRL) to the plasma membrane. Acts as a receptor for adrenomedullin (AM) together with CALCRL. [UniProtKB/Swiss-Prot Function]



[View online »](#)