

## Product datasheet for **MR201812L4V**

### 2010315L10Rik (BC038519) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	2010315L10Rik (BC038519) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	2010315L10Rik
Synonyms:	2010315L10Rik; 5730403H22Rik; AV002165; D12; Ed2; p31; Q-snare
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC038519
ORF Size:	570 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201812).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC038519</a> , <a href="#">AAH38519</a>
RefSeq Size:	933 bp
RefSeq ORF:	572 bp
Locus ID:	67023
Cytogenetics:	8 B3.3
Gene Summary:	SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.[UniProtKB/Swiss-Prot Function]



[View online »](#)