

Product datasheet for **MR201823L4V**

Cbx5 (NM_001110216) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cbx5 (NM_001110216) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cbx5
Synonyms:	2610029O15Rik; C75991; HP1; Hp1a; Hp1alpha
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001110216
ORF Size:	576 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201823).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001110216.1 , NP_001103686.1
RefSeq Size:	8902 bp
RefSeq ORF:	576 bp
Locus ID:	12419
UniProt ID:	Q61686
Cytogenetics:	15 F3



[View online »](#)

Gene Summary:

Component of heterochromatin that recognizes and binds histone H3 tails methylated at 'Lys-9' (H3K9me), leading to epigenetic repression. In contrast, it is excluded from chromatin when 'Tyr-41' of histone H3 is phosphorylated (H3Y41ph). Can interact with lamin-B receptor (LBR). This interaction can contribute to the association of the heterochromatin with the inner nuclear membrane. Involved in the formation of functional kinetochore through interaction with MIS12 complex proteins (By similarity).[UniProtKB/Swiss-Prot Function]