

Product datasheet for **MR201825L4V**

Adat2 (NM_025748) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adat2 (NM_025748) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adat2
Synonyms:	4933426M09Rik; Deadc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025748
ORF Size:	576 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201825).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025748.3
RefSeq Size:	1453 bp
RefSeq ORF:	576 bp
Locus ID:	66757
UniProt ID:	Q6P6J0
Cytogenetics:	10 A2
Gene Summary:	Probably participates in deamination of adenosine-34 to inosine in many tRNAs. [UniProtKB/Swiss-Prot Function]



[View online »](#)