

Product datasheet for **MR201849L3V**

Nkiras1 (NM_023526) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Nkiras1 (NM_023526) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Nkiras1 |
| Synonyms: | 2400004O09Rik; kappaB-Ras1; KBRAS1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_023526 |
| ORF Size: | 579 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR201849). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_023526.3 , NP_076015.2 |
| RefSeq Size: | 4592 bp |
| RefSeq ORF: | 579 bp |
| Locus ID: | 69721 |
| UniProt ID: | Q8CEC5 |
| Cytogenetics: | 14 A1 |



[View online »](#)

Gene Summary:

Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and mediating cytoplasmic retention of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB (By similarity).[UniProtKB/Swiss-Prot Function]