

## Product datasheet for MR201853

### Sec11c (NM\_025468) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sec11c (NM\_025468) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Sec11c  
**Synonyms:** 1810029G24Rik; Sec1113  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201853 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGTTCGTGCGGGCGCCGTGGGACGCATCTCCCACGTCCAGCCTGGACATCTTTGGGGACCTGAGGA  
 AGATGAACAAACGACAGCTCTACTACCAGGTTTTAAACTTTGCCATGATCGTGTCTTCTGCGCTCATGAT  
 CTGGAAAGGCCTGATTGTTCTCACGGGCAGCGAGAGTCCCATCGTGGTGGTACTCAGTGGCAGTATGGAG  
 CCGGCCTCCACAGAGGAGATCTGCTGTCTCACGAATTTCCGGGAGGACCCCATCAGAGCTGGTGAAA  
 TAGTTGTTTTTAAGGTTGAAGGAAGAGACATTCGATAGTTCACAGAGTAATCAAGTTTCATGAAAAAGA  
 TAATGGTGACATCAAGTTTCTGACTAAAGGAGATAATAATGAAGTCGATGATAGAGGCTTGACAAAGAA  
 GGCCAGAACTGGCTGGAGAAGAAGGACGTGGTAGGAAGAGCTCGAGGGTCTTACCTTATGTTGGCATGG  
 TCACCATCATAATGAATGACTATCCAAAATTCAAGTATGCTCTTTGGCTGTGATGGGTGCATATGTGTT  
 ACTAAAACGTGAATCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR201853 protein sequence  
 Red=Cloning site Green=Tags(s)

MVRAGAVGTHLPTSSLDIFGDLRKNKRQLYYQVLNFMIVSSALMIWKGLIVLTGSESPIVVVLGSME  
 PAFHRGDLLFLTNFREDPIRAGEIVVFKVEGRDIPVHRVIKVEKDNKFLTKGDNNEVDDRGLYKE  
 GQNWLEKKDVVGRARGFLPYVGMVTIIMNDYPKFKYALLAVMGAYVLLKRES

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV



Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM\_025468

ORF Size: 579 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_025468.2](#)

RefSeq Size: 770 bp

RefSeq ORF: 579 bp

Locus ID: 66286


UniProt ID: [Q9D8V7](#)

**Cytogenetics:** 18 E1

**MW:** 21.7 kDa

**Gene Summary:** Component of the microsomal signal peptidase complex which removes signal peptides from nascent proteins as they are translocated into the lumen of the endoplasmic reticulum.  
[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR201853