

## Product datasheet for **MR201853L3V**

### Sec11c (NM\_025468) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sec11c (NM_025468) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sec11c
Synonyms:	1810029G24Rik; Sec11I3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025468
ORF Size:	579 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201853).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025468.1</a>
RefSeq Size:	770 bp
RefSeq ORF:	579 bp
Locus ID:	66286
UniProt ID:	<a href="#">Q9D8V7</a>
Cytogenetics:	18 E1
Gene Summary:	Component of the microsomal signal peptidase complex which removes signal peptides from nascent proteins as they are translocated into the lumen of the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]



[View online »](#)