

Product datasheet for **MR201864A1V**

Mouse *Csrp2* (NM_007792) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse *Csrp2* (NM_007792) AAV Particle
Tag: Myc-DDK
Symbol: *Csrp2*
Synonyms: AW551867; Crp2; SmLim
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201864 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCTGTCTGGGGCGGTGGAATAAGTGCAGGGCCTGCGGGAGAACCGTGTACCACGCGGAAGAGGTGC
 AGTGCGATGGGCGGAGCTTCCATCGCTGCTGCTTCCTGTGCATGGTTGCAGGAAAAATTTAGACAGCAC
 AACAGTGGCGATTTCATGATGAAGAGATCTACTGCAATCCTGCTACGGAAAGAAGTATGGACCAAAAGGC
 TATGGTTATGGCCAGGGCGCTGGCAGCTCAACATGGACCGCGGTGAGAGACTGGGCATCAAGCCAGAGA
 GTGCTCAACCTCACAGGCCTACGACAAATCCAAACACTTCTAAATTTGCCAGAAATATGGAGGAGCTGA
 GAAGTGTTCAGGTGTGGGATTCCGTGTATGCTGCGGAGAAGATCATTGGAGCTGGGAAGCCCTGGCAC
 AAAAACTGTTCCGGTGTGCCAAGTGTGGGAAGAGTCTGGAGTCTACAACTCTGACTGAGAAAGAAGGCG
 AAATCTACTGTAAAGGTGCTACGCAAAGAACTTTGGGCCAAGGGATTTGGCTATGGTCAAGGGGCAGG
 GGCCCTTGTTTCATGCTCAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201864 protein sequence
 Red=Cloning site Green=Tags(s)

MPVWGGGNKCGACGRVYHAEVQCDGRSFHRCFLCMVCRKNLDSTTVAIHDEEYCKSCYKKGYPKG
 YGYGQGAGTLNMDRGERLGIKPESAQPHRPTTNPNTSKFAQKYGGAEKCSRCGDSVYAAEKIIGAGKPWH
 KNCFRCAKCGKSLESTLTKEGEYCKGKYAKNFGPKGFGYQGAGALVHAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_007792
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_007792.2
RefSeq Size:	889 bp
RefSeq ORF:	582 bp
Locus ID:	13008
UniProt ID:	P97314
Cytogenetics:	10 D1
MW:	20.9 kDa