

Product datasheet for MR201873

Csrp1 (NM_007791) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Csrp1 (NM_007791) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Csrp1
Synonyms: AA408841; AA959891; AW545626; CRP1; Csrp
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201873 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCAACTGGGGAGGAGGCAAGAAATGTGGGTATGCCAGAAGACGGTCTACTTTGCTGAGGAGGTCC
 AGTGCGAGGGCAACAGCTTCCATAAATCCTGCTTCTGTGTATGGTCTGCAAGAAAAATCTGGACAGCAC
 CACTGTGGCAGTGCATGGAGAAGAGATCTATTGCAAGTCATGTTACGGCAAGAAGTACGGGCCAAAAGGC
 TACGGCTACGGGCAGGGCGCAGGCACGCTGAGCACAGACAAGGGGAGTCTCTGGGCATCAAGCATGAGG
 AAGCCCTGGACACAGGCCACCACCAACCCCAATGCATCCAAGTTTGGCCAGAAGATTGGCGGCTCTGA
 GCGCTGTCCCGCTGTAGCCAGGCGGTCTATGCTGCGGAGAAGGTGATCGGTGCAGGGAAGTCTGGCAT
 AAGTCTGCTTCCGATGTGCAAGTGTGGCAAAGGCCCTTGAGTCGACCACCCTGGCAGACAAGGATGGTG
 AGATCTACTGTAAAGGATGCTATGCCAAAACTTTGGGCCAAAGGTTTTGGCTTTGGACAGGGAGCTGG
 AGCCTTGGTTCACTCAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201873 protein sequence
 Red=Cloning site Green=Tags(s)

MPNWGGGKKKGCVCQKTVYFAEEVQCEGNSFHKSCFLCMVCKKNLDSTTVAVHGEEIYCKSCYKGYGPKG
 YGYGQAGTLSTDKGESLGIKHEEAPGHRPTTNPNAKFAQKIGGSERCPRCSQAVYAAEKVIGAGKSWH
 KSCFRCAKCGKLESTTLADKDEIYCKGKYAKNFGPKGFVGGQAGALVHSE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_007791

ORF Size: 582 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007791.1](#)

RefSeq Size: 1770 bp

RefSeq ORF: 582 bp

Locus ID: 13007


UniProt ID: [P97315](#)

Cytogenetics: 1 E4

MW: 20.6 kDa

Gene Summary: Could play a role in neuronal development.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201873