

Product datasheet for **MR201878L3V**

Ak1 (BC014802) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ak1 (BC014802) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ak1
Synonyms:	Ak-1; B430205N08Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC014802
ORF Size:	582 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201878).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC014802 , AAH14802
RefSeq Size:	2059 bp
RefSeq ORF:	584 bp
Locus ID:	11636
Cytogenetics:	2 22.09 cM



[View online »](#)

Gene Summary:

Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. Also possesses broad nucleoside diphosphate kinase activity. Plays an important role in cellular energy homeostasis and in adenine nucleotide metabolism (By similarity). May provide a mechanism to buffer the adenylate energy charge for sperm motility.
[UniProtKB/Swiss-Prot Function]