

Product datasheet for **MR201892A1V**

Mouse Rnf41 (BC019415) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Rnf41 (BC019415) AAV Particle
Tag: Myc-DDK
Symbol: Rnf41
Synonyms: 4933415P08Rik, FLRF
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201892 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCCAAAGATGAAGTCCCAAACCACAATTGCATTAAGCACCTGCGCTCCGTGGTCCAGCAGCAGCAGT
CGCGCATCGCAGAGCTGGAGAAGACCAGGGCTGAACACAAGCACCAGCTGGCAGAGCAGAAGCGAGACAT
TCAGCTGCTGAAGGCGTATATGCGAGCCATCCGCACTGTCAACCCCAACCTTCAGAACCTGGAGGAGACA
ATCGAATAACAAGAGATCCTCGAGTGGGTGAATCCCTGCAGCCGGCAAGGGTGACCCGCTGGGGGGCA
TGATCTCCACTCCTGATGCTGTGCTCCAGGCTGTCATCAAGCGCTCCCTCGTGGAAAGTGGCTGCCCGGC
CTCCATCGTCAACGAGCTGATTGAAAATGCCCATGAACGCACTTGGCCCCAGGGTCTGGCCACACTAGAG
ACAAGACAGATGAACCGCGCTACTATGAGAACTACGTGGCCAAGCGCATCCCTGGCAAGCAGGCTGTAG
TGGTGATGGCCTGTGAGAACCAGCACATGGGGACGACATGGTGCAGGAGCCAGGGCTCGTCATGATATT
TGCGCATGGTGTGGAGGAGATA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201892 protein sequence
Red=Cloning site Green=Tags(s)

MPKDELPHNHCILKLSRVVQQQSRIAELEKTRAHKKHLAEQKRDIQLLKAYMRAIRSVNPNLQNLLEET
IEYNEILEWVNSLQPARVTRWGMISTPDAVLQAVIKRSLVESGCPASIVNELIENAHERSWPQGLATLE
TRQMNRYYENYVAKRIPGKQAVVVMACENQHMGDDMVQEPGLVMIFAHGVVEEI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	BC019415
ORF Size:	582 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC019415 , AAH19415
RefSeq Size:	1739 bp
RefSeq ORF:	584 bp
Locus ID:	67588
Cytogenetics:	10 76.55 cM
MW:	22.3 kDa