

## Product datasheet for MR201893L3

### Rnf4 (NM\_011278) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf4 (NM_011278) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Rnf4
Synonyms:	AU018689; Gtrge; Gtrgeo8; SNURF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201893).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_011278
ORF Size:	582 bp



[View online >](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_011278.4</a>
<b>RefSeq Size:</b>	2936 bp
<b>RefSeq ORF:</b>	585 bp
<b>Locus ID:</b>	19822
<b>UniProt ID:</b>	<a href="#">Q9QZS2</a>
<b>Cytogenetics:</b>	5 B2
<b>Gene Summary:</b>	This gene encodes a small nuclear RING finger protein that mediates ubiquitylation of polysumoylated proteins. Deficiency of the encoded protein in mice leads to embryonic lethality and global DNA hypermethylation. A similar protein in humans is required for arsenic-induced degradation of promyelocytic leukemia protein in acute promyelocytic leukemia. Alternative splicing of this gene results in multiple transcript variants. A pseudogene for this gene has been identified on chromosome 10. [provided by RefSeq, Jan 2015]

Product images:



Circular map for MR201893L3