

Product datasheet for **MR201903L4V**

Lman2l (BC027226) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lman2l (BC027226) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lman2l
Synonyms:	VIPL, VIP36-like
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC027226
ORF Size:	585 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201903).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC027226.1
RefSeq Size:	2199 bp
RefSeq ORF:	587 bp
Locus ID:	214895
Cytogenetics:	1 B
Gene Summary:	May be involved in the regulation of export from the endoplasmic reticulum of a subset of glycoproteins. May function as a regulator of ERGIC-53 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)