

Product datasheet for **MR201913A1V**

Mouse Bid (NM_007544) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Bid (NM_007544) AAV Particle
Tag: Myc-DDK
Symbol: Bid
Synonyms: 2700049M22Rik; AI875481; AU022477
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201913 representing NM_007544
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGACTCTGAGGTCAGCAACGGTTCGGCCTGGGGGCCGAGCACATCACAGACCTGCTGGTGTTCCGGCT
 TTCTCCAAAGCTCTGGCTGTACTCGCCAAGAGCTGGAGGTGCTGGGTCGGGAAGTGCCTGTGCAAGCTTA
 CTGGGAGGCAGACCTCGAAGACGAGCTGCAGACAGATGGCAGCCAGGCCAGCCGCTCCTTCAACCAAGGA
 AGAATAGAGCCAGATTCTGAAAGTCAGGAAGAAATCATCCACAACATTGCCAGACATCTCGCCCAATAG
 GCGATGAGATGGACCACAACATCCAGCCACACTGGTGAGACAGCTAGCCGCACAGTTCATGAATGGCAG
 CCTGTCCGAGGAAGACAAAAGGAAGTGCCTGGCCAAAGCCCTTGATGAGGTGAAGACAGCCTTCCCCAGA
 GACATGGAGAACGACAAGGCCATGCTGATAATGACCATGCTGTTGGCCAAAAAAGTGGCCAGTCACGCAC
 CATCTTTGCTCCGTGATGTCTTCCACACGACTGTCAACTTTATTAACCAGAACCTATTCTCCTATGTGAG
 GAACTTGGTTAGAAACGAGATGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201913 representing NM_007544
 Red=Cloning site Green=Tags(s)

MDSEVSNGLGAEHITDLLVFGFLQSSGCTRQELEVLGRELVPQAYWEADLEDELQTDGSQASRSFNQG
 RIEPDSESQEEIIHNIARHLAQIGDEMDHNIQPTLVRQLAAQFMNGSLSEEDKRNCLAKALDEVKTAFFR
 DMENDKAMLIMTMLLAKKVASHAPSLLRDVFHTTVNFINQNLFSYVRNLVRNEMD

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_007544
ORF Size:	585 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_007544.3
RefSeq Size:	1906 bp
RefSeq ORF:	588 bp
Locus ID:	12122
UniProt ID:	P70444
Cytogenetics:	6 57.02 cM
MW:	22 kDa