

Product datasheet for MR201913L4V

Bid (NM_007544) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bid (NM_007544) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bid
Synonyms:	2700049M22Rik; AI875481; AU022477
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007544
ORF Size:	588 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201913).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007544.3
RefSeq Size:	1906 bp
RefSeq ORF:	588 bp
Locus ID:	12122
UniProt ID:	P70444
Cytogenetics:	6 57.02 cM
Gene Summary:	Induces caspases and apoptosis. Counters the protective effect of Bcl-2. The major proteolytic product p15 BID allows the release of cytochrome c.[UniProtKB/Swiss-Prot Function]



[View online »](#)