

Product datasheet for **MR201916L3V**

Nmrk1 (NM_145497) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nmrk1 (NM_145497) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nmrk1
Synonyms:	AA408380; D630020N23Rik; Nrk1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145497
ORF Size:	588 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201916).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145497.2 , NP_663472.1
RefSeq Size:	1637 bp
RefSeq ORF:	588 bp
Locus ID:	225994
UniProt ID:	Q91W63
Cytogenetics:	19 B
Gene Summary:	Catalyzes the phosphorylation of nicotinamide riboside (NR) and nicotinic acid riboside (NaR) to form nicotinamide mononucleotide (NMN) and nicotinic acid mononucleotide (NaMN). [UniProtKB/Swiss-Prot Function]



[View online »](#)