

## Product datasheet for **MR201922L4V**

### **Bcl2l11 (NM\_207680) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bcl2l11 (NM_207680) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bcl2l11
Synonyms:	1500006F24Rik; bcl2-L-11; Bim; Bod
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207680
ORF Size:	591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201922).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207680.2</a> , <a href="#">NP_997563.1</a>
RefSeq Size:	5040 bp
RefSeq ORF:	591 bp
Locus ID:	12125
UniProt ID:	<a href="#">O54918</a>
Cytogenetics:	2 F1
Gene Summary:	Induces apoptosis and anoikis. The isoforms vary in cytotoxicity with isoform BimS being the most potent and isoform BimEL being the least potent.[UniProtKB/Swiss-Prot Function]



[View online »](#)