

## Product datasheet for **MR201949L3V**

### Stom (BC003789) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Stom (BC003789) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Stom
Synonyms:	Epb7.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC003789
ORF Size:	591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201949).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC003789</a> , <a href="#">AAH03789</a>
RefSeq Size:	2191 bp
RefSeq ORF:	593 bp
Locus ID:	13830
Cytogenetics:	2 23.53 cM
Gene Summary:	Regulates ion channel activity and transmembrane ion transport. Regulates ASIC2 and ASIC3 channel activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)